




Complete Gas Analyser



Scion 8500-GC

- Simultaneous analysis of He, H₂, O₂, N₂, CH₄, CO, CO₂, C₂H₄, C₂H₆, C₂H₂, H₂S from 50ppm to 100% and C₃ to C₁₂ from 1ppm to 100%
- Less than 4 minutes up to C₆ and < 10 minutes to C₁₂
- Any sample pressure –atmospheric to 1200psi
- Moisture is no problem – columns are backflushed
- Single combined report for all analytes
- Optional software-controlled stream selection
- Optional oxygenates, alkenes, alkynes and aromatics



Custom Analytical Solutions for Gas Analysis

Contact sales-can@scioninstruments.com for more information